



ELECTROSURGICAL UNITS

KENTAMED 1E



APPLICATION: MINOR and MEDIUM SURGERY

USE: In medical and dental offices and clinics

DESCRIPTION:

- ❖ intelligent, functionally complete, microprocessor controlled unit
- ❖ 100 W maximal output power
- ❖ 5+1 operating modes:
 - smooth cutting
 - cutting with increased degree of coagulation
 - contact coagulation
 - spray coagulation
 - bipolar and micro bipolar coagulation
- ❖ independent power setting for each mode – digital LED
- ❖ digital memory for last power setting
- ❖ overheating (OH) and power failure (PF) sensors
- ❖ status monitor on the front panel
- ❖ automatic and immediate output power switch off if controlled parameters are out of norm
- ❖ sound alarm if controlled parameters are out of norm
- ❖ 3 year warranty

KENTAMED 1E

Technical Characteristics		
Mode	Maximal Output Power	Crest Factor
Pure CUT	100 W / 500 Ω	1.6
Blend CUT	100 W / 500 Ω	2.1
Contact COAG	100 W / 300 Ω	3.0
Spray COAG	20 W / 5000 Ω	7.0
Bipolar COAG	80 W / 100 Ω	1.5
Micro Bipolar COAG	35 W / 100 Ω	1.5
Operating Frequency [kHz]	460	
Modulation Frequency [kHz]	33	
Patient Safety	NEMSY (Neutral Electrode Monitoring System)	
Activation	Twin Footswitch and Hand-switch Pencil	
Duty Cycle	10/30	
Safety Class/MDD Class	I/IIb	
Applied Parts Type	BF	
Dimensions [mm]	256x200x85	
Weight [kg]	5	

STANDARD SET OF ACCESSORIES:

- ✓ 1 reusable finger-switch handle with 3 m cable
- ✓ 7 monopolar electrodes
- ✓ 1 reusable neutral electrode with 3 m cable
- ✓ 1 bipolar forceps
- ✓ 1 bipolar cable 3 m
- ✓ 1 twin footswitch
- ✓ user manual

OTHER ACCESSORIES – available upon request

Manufacturer:

KENTAMED Ltd.
 41, Nestor Abadjiev St.
 POB 203, Plovdiv 4023, Bulgaria
 tel./fax: +359 32 622 934
 e-mail: sales@kentamed.com
www.kentamed.com



Distributor:

