



THOR LAPAROSCOPIC ELECTRODES®

The new name for quality and precision
in electrosurgical instruments



THOR  **ELECTRODES®**

PTFE Coated Laparoscopic Electrodes



IGS 0017 THOR PTFE Laparoscopic Blade 33cm
IGS 0017S THOR PTFE Laparoscopic Blade 33cm Split Stem



IGS 0018 THOR PTFE Laparoscopic Spatula 33cm
IGS 0018S THOR PTFE Laparoscopic Spatula 33cm Split Stem



IGS 0019 THOR PTFE Laparoscopic Curved Blade 33cm
IGS 0019S THOR PTFE Laparoscopic Curved Blade 33cm Split Stem



IGS 0020 THOR PTFE Laparoscopic L-Hook 33cm
IGS 0020S THOR PTFE Laparoscopic L-Hook 33cm Split Stem
IGS 0020X THOR PTFE Laparoscopic L-Hook 45cm
IGS 0020XS THOR PTFE Laparoscopic L-Hook 45cm Split Stem



IGS 0024 THOR PTFE Laparoscopic Ball 33cm
IGS 0024S THOR PTFE Laparoscopic Ball 33cm Split Stem



IGS 0027 THOR PTFE Laparoscopic Thin Curved 33cm
IGS 0027S THOR PTFE Laparoscopic Thin Curved 33cm Split Stem



IGS 0100 THOR PTFE Laparoscopic Wire L-Hook 33cm
IGS 0100S THOR PTFE Laparoscopic Wire L-Hook 33cm Split Stem
IGS 0100X THOR PTFE Laparoscopic Wire L-Hook 45cm
IGS 0100XS THOR PTFE Laparoscopic Wire L-Hook 45cm Split Stem



IGS 0105 THOR PTFE Laparoscopic Wire J-Hook 33cm
IGS 0105S THOR PTFE Laparoscopic Wire J-Hook 33cm Split Stem

FEATURES & BENEFITS

✓ THOR Laparoscopic Electrodes are coated in high-grade PTFE, making surgery faster and less disrupted.

✓ Reintroduction during surgery is minimised whilst lower power settings make for safer surgery.

✓ Disposable electrodes provide a new fully insulated shaft for each case, preventing stray current and unintended burns.



IGS 0021 THOR PTFE Laparoscopic J-Hook 33cm
IGS 0021S THOR PTFE Laparoscopic J-Hook 33cm Split Stem
IGS 0021X THOR PTFE Laparoscopic J-Hook 45cm
IGS 0021XS THOR PTFE Laparoscopic J-Hook 45cm Split Stem



IGS 0022 THOR PTFE Laparoscopic Needle 33cm
IGS 0022S THOR PTFE Laparoscopic Needle 33cm Split Stem



THOR ELECTRODES®



IG Surgical Generator



Standard End - Compatible with hand switching



Split Stem - Compatible with foot switching



GET IN TOUCH!

THOR ELECTRODES®



+ 44 (0) 1242 262680



enquiries@igsurgical.co.uk



www.igsurgical.co.uk



IG Surgical Office:

Inter Global Surgical LLP
Suite 100
20 Winchcombe Street
Cheltenham
GL52 2LY

Additional Items



Our smoke evacuation solution works for both open and closed procedures and with any ESU.



The Auckland Patient Return Electrodes improve safety, increase simplicity and save money. Full two year warranty.

Please contact us for more information or to arrange a free trial.