

MicroTip™ Plus

Ultrasonic Surgical Tip



Because we are committed to limit uncertainty, the Integra MicroTip™ Plus ultrasonic surgical tip provides surgeons with a safe and comfortable solution for transnasal procedures.

Product Description


- Improved access for transnasal and other procedures requiring extended length.
- Small diameter for improved visualization and access to the surgical site.
- Functional with TissueSelect™ technology.
- Tip changeable in sterile field.
- Pre-aspiration holes cool the tip and reduce misting.
- Available on the small and powerful 36 kHz handpiece for improved tactile feedback.



Compatible with CUSA EXcel+
Ultrasonic Tissue Ablation System

MicroTip™ Plus

Ultrasonic Surgical Tip



Product Specifications

| | |
|--|---------------------|
| Catalog number | C4615S |
| Tip Length | 19.27 cm (7.59 in.) |
| Tip ID | 1.57 mm (0.06 in.) |
| Packaging: 4 single use sterile tips per box | |



CUSA EXcel® + Range of Tips

PrecisionTip™ Series

Smallest diameter for access in restricted areas



MicroTip™ Series

Better access for delicate procedures



Standard Tip Series

Larger diameter for rapid tumor debulking



ShearTip™

Specific design for fibrous and tough tissue



CUSA® SaberTip™

Selective bone removal around critical structures



MicroTip™ Plus

Improved access for transnasal procedures



Integra LifeSciences Services (France) SAS
Sales & Marketing EMEA
Immeuble Séquoia 2 • 97 allée Alexandre Borodine
Parc technologique de la Porte des Alpes
69800 Saint Priest • FRANCE
T +33 (0)4 37 47 59 00 • fax +33 (0)4 37 47 59 99
emea.info@integralife.com • integralife.com

Customer Service

International: +33 (0)4 37 47 59 50 • +33 (0)4 37 47 59 25 (Fax) • csema@integralife.com
France: +33 (0)4 37 47 59 10 • +33 (0)4 37 47 59 29 (Fax) • custsvfrance@integralife.com
United Kingdom: +44 (0)1 264 345 780 • +44 (0)1 264 363 782 (Fax) • custsvuk@integralife.com
Germany: +49 (0) 2102 5535 6200 • +49 (0)2 102 5536 636 (Fax) • custsvgermany@integralife.com
Benelux: +32 (0)2 257 4130 • +32 (0)2 253 2466 (Fax) • custsvbenelux@integralife.com
Switzerland: +41 (0)2 27 21 23 30 • +41 (0)2 27 21 23 99 (Fax) • custsvcsuisse@integralife.com

Integra LifeSciences (Ireland) Limited
IDA Business and Technology Park
Sragh
Tullamore • County Offaly • Ireland

©2011 Integra LifeSciences Corporation. All rights reserved. ILS 04-01-004-02-11
PRODUCTS FOR SALE IN EUROPE, MIDDLE-EAST AND AFRICA ONLY

Distributed by

INTEGRA
LIMIT UNCERTAINTY

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region. • Always refer to the appropriate instructions for use for complete clinical instructions. • Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality. • WARNING: Applicable laws restrict these products to sale by or on the order of a physician. • CUSA EXcel, Microtip, Precision Tip, ShearTip, SaberTip, MicroTip Plus, Integra and the Integra logo are trademarks or registered trademarks of Integra LifeSciences Corporation or its subsidiaries.



CE
0086

IMPRIM'VERT™