329522

MEDICATION TROLLEY, 45+45 CM DRAWERS, WITH UPPER ACCESSORY-HOLDER STRUCTURE - LINKAR RANGE

MD10 Rev. 06 Data: 11/09/24 - indicative image



The structural body of the trolley is made by assembling composite materials, each one of which is chosen to optimise the functionality, reliability and sanitising of the product. Defined by an attentive ergonomic study, it guarantees maximum stability in every condition of use and offers excellent manoeuvrability, even when fully loaded. Self-bearing structure, made of epoxy powder coated steel sheet with antibacterial protection that ensures a higher degree of hygiene and sanitation. The trolley, characterised by rounded lines, without corners and sharp edges, allows optimum sanitation.
The drawers are made of epoxy powder-coated sheet steel, with antibacterial additive, and an injection-moulded plastic front with rounded edges and integrated coloured handles, complete with label holders. Easy-grip ergonomic handles. The drawers are fitted on total extraction steel guides, with spring damper closure and lever-operated unhooking system. The drawers, of each compartment, are closed by
means of a centralized locking system with folding-anti-trauma key, located on the front.
The right side is provided with 2 metal accessory-holder bars.
The left side is equipped with an additional pull-out worktop and with 2 metal accessory-holder bars.
The top is made of thermoformed ABS plastic, impact- and scratch-resistant, easy to wash and disinfect. The top is complete with a raised fluid retaining front edge and a built-in raiser on three sides which can be used as a push handle. The internal shaping of the built-in raiser provides an easy grip and enhances the manoeuvrability of the trolley.
The base, made of impact- and scratch-resistant thermoformed ABS plastic, is easy to wash and disinfect. Complete with plastic bumper wheels, it is equipped with 4 swivel wheels with a 125 mm diameter, of which 2 at the front with a foot-operated brake, assembled onto a steel frame.
The trolley is supplied as standard with:
3+3 drawers width 450 mm and h 150 mm, blue handle
1+1 drawer width 450 mm and h 225 mm, green handle
1 upper accessory-holder structure with two bars, the lower adjustable
1 kit of 9 flap-down drawers (5+4)
2 large stainless steel shelves with 1 plastic divider
1 double hook for I.V. for upper structure
1 support with stainless steel kidney dish
1 disinfection basin with support for upper structure
1 universal support
1 yellow container for small waste
Worktop h 101 cm
Overall dimensions: cm 117x62,5x172 H
Safe working load: 100 kg
MEDICAL DEVICE, CLASS I, COMPLIANT WITH REGULATION (EU) 2017/745
MEDICAL DEVICE MANUFACTURED IN COMPLIANCE WITH ISO 9001, ISO 13485 AND ISO 14001 CERTIFICATION

It is not authorized coping, however in part, or using of the present drawing.The Company reserves the right to introduce all technical and/or aesthetical modifications deemed necessary to their products without prior notice.

Malvestio S.P.A. 35010 Via Marconi 12/D, Loc. Murelle, 35010 Villanova di Camposampiero (PD) Italia
www.malvestio.it - Tel. 049/9299511 ? Fax 049/9299500 ? Mail:info@malvestio.it ? cap. soc. 2.970.000,00 ? i.v. ? c.c.i.a.a. Padova 119204 ? reg.imprese, c.f. e p.iva 00197370281- reg. AEE IT22030000013732

Only products that are identified as such are FSC® certified