32960

BASIC TROLLEY WITH 9 GUIDES, 60 CM DRAWERS, 2 ACCESSORY HOLDER SIDE BARS RIGHT SIDE, PULL-OUT WORKTOP AND 3 FLAP-DOWN DRAWERS LEFT SIDE - LINKAR RANGE

MD10 Rev. 01 data: 10/09/24 - indicative image



The trolley is made by assembling composite materials, each one of which is chosen to optimise the functionality, reliability and sanitising of the product. Defined by an attentive ergonomic study, it guarantees maximum stability in all condition of use and offers excellent manoeuvrability, even when fully loaded.
Self-bearing structure made of epoxy powder coated with antibacterial protection. It provides durable and effective against harmful bacteria, moulds and fungi, reducing the risk of cross-contamination and complementing existing hygiene protocols.
The trolley, characterised by rounded lines, without corners and sharp edges, allows optimum sanitation.
Inner sides can be equipped with drawers (accessories), 60 cm wide, complete with total extraction steel guides, available in 3 heights (7,5 - 15 - 22,5 cm), and with removable trays (accessories), 60 cm wide, for medication bins.
Drawers and trays are closed by means of a centralized locking system with folding-anti-trauma key, located on the front edge.
The right side is equipped with 2 metal accessory-holder bars.
The left side is equipped with an additional pull-out worktop and 3 flap-down drawers, in translucent plastic, with a 40° opening angle, which can be removed for cleaning and sanitisation.
The top is made of thermoformed ABS plastic, impact- and scratch-resistant, easy to wash and disinfect. The top is complete with a raised fluid retaining front edge and a built-in raiser on three sides which can be used as a push handle. The internal shaping of the built-in raiser provides an easy grip and enhances the manoeuvrability of the trolley. The trolley can be equipped with an upper accessory-holder structure, either h 30 cm or h 60 cm, fixed on the rear corners of the top, or with other supports.
The base made of thermoformed ABS plastic, impact- and scratch-resistant, is easy to clean and sanitise. Complete with plastic bumper wheels, it is equipped with 4 swivel wheels, diameter 125 mm, 2 of which, at the front, with a foot-operated brake, assembled onto a steel frame.
Worktop h 101 cm
Overall dimensions: cm 91x62,5x105,5 H
Safe working load: 100 kg
MEDICAL DEVICE, CLASS I, COMPLIANT WITH REGULATION (EU) 2017/745
MEDICAL DEVICE MANUFACTURED IN COMPLIANCE WITH ISO 9001, ISO 13485 AND ISO 14001 CERTIFICATION.

MALVESTIO®
Made in Italy

It is not authorized coping, however in part, or using of the present drawing.The Company reserves the right to introduce all technical and/or aesthetical modifications deemed necessary to their products without prior notice.

Malvestio S.P.A. 35010 Via Marconi 12/D, Loc. Murelle, 35010 Villanova di Camposampiero (PD) Italia
www.malvestio.it - Tel. 049/9299511 ? Fax 049/9299500 ? Mail:info@malvestio.it ? cap. soc. 2.970.000,00 ? i.v. ? c.c.i.a.a. Padova 119204 ? reg.imprese, c.f. e p.iva 00197370281- reg. AEE IT22030000013732

Only products that are identified as such are FSC® certified